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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

## INDICATIONS

- Grade II/III ankle sprain.
- Post ligament, soft tissue and tendon surgery.
- Stable fracture of the fibula, talus, calcaneus and medial malleolus.

## PERFORMANCE

The devices provide immobilization/motion restriction of the ankle.

## CONTRAINDICATION

- Paresthesia
- Severe varicose vein
- Skin disorders or injuries in the relevant area
- Impaired lymphatic drainage

## WARNINGS AND PRECAUTIONS OF USE

- The orthosis is a single patient use device. Do not use for multiple patients.
- Do not wear for extended periods without medical follow-up.
- Do not remove the device without medical advice.
- The orthosis must not be used in direct contact with damaged skin.
- If skin irritations of any kind develop, consult medical professional immediately.
- When adjusting the strap, tighten it so that it is firm but comfortable.
- If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend

### A. Donning the device :

1. Unfasten all the straps with hook-and-loop and the hook-and-loop located on the top of the liner. Remove the front panel to pull the textile liner out the front of the walker.
2. Place the foot in the textile lining and wrap the foot and lower leg well (schema 1).
3. Place the wrapped foot in the walking boot. The heel should press against the back of the walker.
4. Place the front panel and make sure the textile lining is well covered (schema 2).
5. Tighten all straps from distal to proximal. Adapt the tightening to obtain a firm but comfortable tightening.

### B. Donning additional foam padding :

1. To account for chafing or spacing, additional foam pads may be placed inside the walker.
2. Secure the pads using the small self-gripping hooks on the side of the leg, between the textile part and the leg, at the level of the plastic side uprights (schema 2). This operation must be carried out when placing the foot in the textile lining (schema 1).

### C. Walker inflation :

To inflate and adjust the air cell pressure, squeeze the bilateral "inflate tip" pump (schema 4 and 5).

### D. Walker deflation :

To deflate the air cells, press the small "deflate tip" button located next to the pump (schema 5).

## STORAGE, MAINTENANCE, DISPOSAL

- Remove straps, pads or liners and close all hook and loop closures before washing. Hand wash liners with mild detergent in cold water and hang dry only. Do not dry clean, bleach and iron.
- The device and packaging must be disposed of in accordance with respective local or national environmental regulation.

## MEASUREMENT AND SIZE CHART

Size of the walker to choose in accordance with shoe size.

Size	2	3	4	5
Shoe size	36-40	40-43	43-46	46-49

## CLASSIFICATION

Class I medical device, non-sterile, no measurement function. This product fulfils the Medical Device Regulation (EU) 2017/745.

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