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## INTRODUCTION SENTENCE

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

## INDICATIONS

- Grade II/III ankle sprain.
- Post ligament, soft tissue and tendon surgery.
- Stable fracture of the fibula, talus, calcaneus and medial malleolus.

## CONTRAINDICATION

- Paresthesia
- Severe varicose vein
- Skin disorders or injuries in the relevant area
- Impaired lymphatic drainage

## WARNINGS AND PRECAUTIONS OF USE

- The orthosis is a single patient use device. Do not use for multiple patients.
- Do not wear for extended periods without medical follow-up.
- Do not use more than 24 hours of continuous.
- The orthosis must not be used in direct contact with a damaged skin.
- If allergic skin rashes occur, please consult your physicians.
- When adjusting the strap, tighten it so that it is firm but comfortable.
- If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## DONNING AND DOFFING


The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend

### A. Donning the walker :

- Unfasten all the straps, remove front and back panels, and open liner. Do not release the straps off the strut (picture 1).
- Place foot and leg inside foam liner and wrap foot and lower leg snugly (picture 2).
- Slip foot inside the walking boot, be sure to place heel close to the heel frame, and make sure the foot rests completely flat down in the brace.
- Close the hook and loop straps from bottom to the top. Start at the foot and working up the leg (picture 3). The straps should be snug enough but not so tight as to restrict circulation or cause discomfort.

### B. Donning additional foam padding :

- To account for chafing or spacing, additional foam pads may be placed inside the boot.
- Put additional foam padding around ankle and foot to fit completely (picture 2).

## STORAGE, MAINTENANCE, DISPOSAL

- Remove straps, pads or liners and close all hook and loop closures before washing. Hand wash liners with mild detergent in cold water and hang dry only. Do not dry clean, bleach and iron.
- The device and packaging must be disposed of in accordance with respective local or national environmental regulation.

## MEASUREMENT AND SIZE CHART

Size of the walker to choose in accordance with shoe size.

Size	1	2	3	4	5
Shoe size	32-36	36-40	40-43	43-46	46-49

## CLASSIFICATION

Class I medical device, non-sterile, no measurement function. This product fulfils the Medical Device Regulation (EU) 2017/745.

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